متن نشانه هایی که برای گزارش آزمایش مایع منی استفاده می شوند (به ترتیب الفبا)

نشانه	شرح
Abnorm	Abnormal:
Brom	Bromelain
EE	Electro-ejaculation
G1	Grade 1; Isolated
G2	Grade 2; Moderate
G3	Grade 3; Large
G4	Grade 4; Gross
Н	Abnormal: High
HGO	Normal; Homogeneous, grey-opalescent
Mec	Mechanically
MIL	Masturbation IN the lab
MOL	Masturbation OUT of the lab
N	Normal
SI	Sexual intercourse using non-toxic condom
VL	Very Low
آخر نمونه	No; Losing last portion
اول نمونه	No; Losing first portion
تبدار	Yes; Febrile
غیر تبدار	Yes; NOT Febrile
كامل	Complete

متن کامنت هایی که در گزارش آزمایش مایع منی استفاده می شوند (به ترتیب الفبا)

عنوان	شرح
Арх	Approximate count: *** spermatozoa in the whole specimen; Motile sp.: ***%; Immotile sp: ***%
C<25; عداد:	Only *** cells were seen in two smears; too few for determination of accurate concentration.
Dif; G: P: L: M:	Differentiation: Immature Germinal cells: ***%; PMNs: ***%; Lymphocytes: ***%; Macrophages: ***%.
n<25; عداد	Only *** Sperms were seen in two preparation; too few for determination of accurate concentration.
n=0	No sperm was seen in two counting preparation; too few for determination of accurate concentration.
No Sperm	No Spermatozoa was seen in the aliquot of the specimen that was assayed.
PH-Vis	Because of high viscosity, the PH value may differ from real value.
Prlg Abs	In prolonged abstinence, the total sperm number dose NOT reflects testicular sperm productivity.
Samp Er:	Because of low concentration, counting precision is low; Sampling error: Approximately ***%.
اثر EE	In Electro-ejaculates, the total sperm number dose NOT reflects testicular sperm productivity.
سانتريفوژ	The specimen was centrifuged to be concentrated; Centrifugation may affect sperm morphology.
شستشو	The specimen was washed to reduce viscosity; Washing may affect sperm morphology.
۶ شمارش	Because the specimen was not homogeneous 3 different duplicate preparations were counted.
میانگین ۶	None of the duplicates were statistically close enough; the count is mean of all the 6 preparations.

هنگام نوشتن جواب بخش های گوناکون آزمایش، به دنبال برخی کامنت ها باید عددهایی مربوط به آن کامنت ها نوشته شود. در متن چنین کامنت هایی سه تا نشانه ی ستاره (***) وجود دارد که هنگام وارد کردن نتیجه در کامپیوتر، آن عددها به جای این ستاره ها گذاشته می شوند.